

## 1. Loading and Cleaning Data:

```
python Copy code  
  
import pandas as pd  
import numpy as np
```

- **Why:** These are necessary libraries for handling data. `pandas` is used for reading and cleaning data, while `numpy` helps manipulate numeric arrays.

```
python Copy code  
  
filenames = ['HR_01.txt', 'HR_02.txt', 'SL_01.txt', 'SL_02.txt', 'SK_01.txt',
```

- **Why:** These represent the EEG data files you are working with. Each file contains raw EEG data for different subjects.

```
python Copy code  
  
def load_eeg_data(filename):  
    data = pd.read_csv(filename, delimiter=',', skiprows=4) # Adjust skiprows  
    numeric_data = data.select_dtypes(include=[np.number])  
    return numeric_data.to_numpy()
```

- **How:** This function reads the data from each file, skips non-relevant rows at the top, and extracts only numeric data columns.
- **Why:** EEG data usually contains both numeric and non-numeric rows (like headers, metadata), which we filter out.

```
python Copy code  
  
eeg_data = {name: load_eeg_data(name) for name in filenames}  
print(eeg_data['HR_01.txt'].shape)
```

- **Why:** Loads and stores EEG data files in a dictionary, and prints the shape to verify data loading.

## 2. Bandpass Filtering:

```
python Copy code  
  
from scipy.signal import butter, filtfilt
```

- **Why:** `butter` designs a Butterworth filter, and `filtfilt` applies it to the data.

```
python Copy code  
  
def butter_bandpass(lowcut, highcut, fs, order=4):  
    nyquist = 0.5 * fs  
    low = lowcut / nyquist  
    high = highcut / nyquist  
    b, a = butter(order, [low, high], btype='band')  
    return b, a
```

- **How:** Creates a bandpass filter to keep frequencies between `lowcut` and `highcut`.
- **Why:** EEG signals have useful frequency ranges (0.5 Hz to 50 Hz), and this filter helps isolate them.

```
python Copy code  
  
def apply_bandpass_filter(data, lowcut, highcut, fs, order=4):  
    b, a = butter_bandpass(lowcut, highcut, fs, order=order)  
    return filtfilt(b, a, data, axis=0)
```

- **Why:** This function applies the bandpass filter to EEG data for noise removal and signal clarity.

```
python Copy code  
  
filtered_eeg_data = {name: apply_bandpass_filter(data, lowcut, highcut, sampling_rate) for name, data in eeg_data.items()}  
print(filtered_eeg_data['HR_01.txt'].shape)
```

- **How:** Filters each EEG data set and stores it.
- **Why:** Ensures only relevant frequencies are retained for later analysis.

### 3. Feature Extraction using FFT:

```
python Copy code  
  
def compute_fft(data, fs):  
    n = len(data)  
    freqs = np.fft.fftfreq(n, 1/fs)  
    fft_values = np.fft.fft(data)  
    fft_magnitude = np.abs(fft_values)  
    return freqs[:n // 2], fft_magnitude[:n // 2]
```

- **How:** Converts time-domain EEG signals into the frequency domain using the Fast Fourier Transform (FFT).
- **Why:** EEG analysis often focuses on frequency components (like alpha, beta waves), which are extracted using FFT.

```
python Copy code  
  
def bandpower(data, freqs, low, high):  
    idx_band = np.logical_and(freqs >= low, freqs <= high)  
    return np.mean(data[idx_band])
```

- **Why:** Calculates the power (strength) of a specific frequency band, allowing us to quantify brainwave activity.

```
python Copy code  
  
bands = {'delta': (0.5, 4), 'theta': (4, 8), 'alpha': (8, 12), 'beta': (12, 30)}  
features = {}  
for name, data in filtered_eeg_data.items():  
    feature_matrix = []  
    for channel in range(data.shape[1]):  
        freqs, fft_magnitude = compute_fft(data[:, channel], sampling_rate)  
        feature_vector = [bandpower(fft_magnitude, freqs, low, high) for band, (low, high) in bands.items()]  
        feature_matrix.append(feature_vector)  
    features[name] = np.array(feature_matrix)
```

- **Why:** Breaks EEG data into different frequency bands (delta, theta, etc.) for each channel and extracts meaningful features from them.

## Feature Extraction with FFT & Power Spectrum:

### 1. Fast Fourier Transform (FFT):

- Converts EEG signals from the time domain to the frequency domain.
- Breaks complex brainwave signals into their component frequencies.
- Allows identification and analysis of brainwave patterns like alpha, beta, delta, theta, and gamma waves.

### 2. Frequency Domain Analysis:

- Enables the study of how much of each frequency is present in the EEG data.
- Reveals specific brainwave activities that correspond to different cognitive or mental states.

### 3. Power Spectrum Analysis:

- Measures the power (strength) of different frequency components.
- Identifies the most prominent brainwave activities by calculating the power in specific frequency bands.

### 4. EEG Frequency Bands:

- **Delta (0.5 – 4 Hz):** Deep sleep, unconsciousness.
- **Theta (4 – 8 Hz):** Light sleep, meditation.
- **Alpha (8 – 12 Hz):** Relaxation, calm focus.
- **Beta (12 – 30 Hz):** Active thinking, concentration.
- **Gamma (30 – 50 Hz):** High-level cognitive processing.

### 5. Bandpower Calculation:

- Computes the average power in specific frequency bands (delta, theta, alpha, beta, gamma).
- Helps quantify brainwave activities and allows comparison between cognitive states.

### 6. Extracted Features:

- Feature matrix includes bandpower values for each EEG channel.
- Useful for machine learning models to classify different mental states.